

DHL current tracking /temperature monitoring devices

- Currently 2 x devices contained in deliveries:
 - MonT2 data logger used to provide a screen for initial identification on whether excursion may have occurred
 - TrackIT is primary monitoring device and temperature reports can be generated remotely by DHL
- Positioned on side of inner carton
- You look in upper left hand corner of the MonT2 screen
 - √ means **temperature in range**
 - **X** means **excursion may have occurred**
- If X displayed vaccines quarantined and NIP Logistics called
- Both devices are returned to courier in shipment container



DHL upgrade to all-in-one tracking/temperature monitoring device

- Deliveries from **28 Nov** will contain TrackIT V3 all-in-one remote tracking/temp monitoring device with LCD screen
- Same positioning on side of inner carton as current devices
- Ambient light on opening delivery box wakes up the screen
- A blue light on the device indicates data is being transmitted to server
- Look in upper left hand corner of screen
 - ✓ means **temperature in range**
 - **X** means **excursion may have occurred**
 - **If screen not working - push button**
- If X is displayed or screen not working after pushing button, quarantine vaccines and call NIP Logistics
- Place device back in shipment container for courier to return to warehouse as usual

